


## NEOTEST™ DOA SALIVA




### INSTRUCTIONS OF USE

#### STEP 1

Saturation indicator  
strip unactivated




OR




At least 3 min.


#### STEP 2



OR



#### STEP 3



Peel off to check the specimen  
while no wicking issue occurred

SALIVA DRUGS OF ABUSE (DOA) COMBINATIONS AVAILABLE


| Product Description                          | Format | Cut-off Value                    | Qualification |
|----------------------------------------------|--------|----------------------------------|---------------|
| 7-Aminonclozepam (ACL) Test                  | Device | 100 ng/mL                        | /             |
| Amphetamine (AMP) Test                       | Device | 50/40 ng/mL                      | CE            |
| Barbiturates (BAR) Test                      | Device | 300/50/30 ng/mL                  | CE            |
| Benzodiazepines (BZO) Test                   | Device | 50/20/10 ng/mL                   | CE            |
| Buprenorphine (BUP) Test                     | Device | 10/5 ng/mL                       | CE            |
| Carisoprodol (SOMA) Test                     | Device | 300 ng/mL                        | /             |
| Cocaine (COC) Test                           | Device | 50/20 ng/mL                      | CE            |
| Codeine (COD) Test                           | Device | 10 ng/mL                         | CE            |
| Cotinine (COT) Test                          | Device | 50/30/10 ng/mL                   | CE            |
| Diphenhydramine (DIP) Test                   | Device | 150/100 ng/mL                    | /             |
| Ecstasy (MDMA) Test                          | Device | 60/50 ng/mL                      | CE            |
| Ethyl Glucuronide (EtG) Test                 | Device | 150/100 ng/mL                    | /             |
| Fentanyl (FEN) Test                          | Device | 10 ng/mL                         | CE            |
| Hydrocodone (HCD) Test                       | Device | 10 ng/mL                         | /             |
| Hydromorphone (HMO) Test                     | Device | 300/150 ng/mL                    | /             |
| Ketamine (KET) Test                          | Device | 100/50 ng/mL                     | CE            |
| Lysergic acid diethylamide (LSD) Test        | Device | 25/10 ng/mL                      | CE            |
| Marijuana (THC) Test                         | Device | 50/40/30/25/15/12/10/5/4/3 ng/mL | CE            |
| Methadone Metabolite (EDDP) Test             | Device | 20 ng/mL                         | CE            |
| Mephedrone (MEP) Test                        | Device | 50 ng/mL                         | /             |
| Methadone (MTD) Test                         | Device | 75/50/30 ng/mL                   | CE            |
| Methamphetamine (MET) Test                   | Device | 50 ng/mL                         | CE            |
| Methaqualone (MQL) Test                      | Device | 150/100 ng/mL                    | CE            |
| Methcathinone (MTC) Test                     | Device | 50 ng/mL                         | /             |
| 3,4-Methylenedioxypropylvalerone (MDPV) Test | Device | 200 ng/mL                        | CE            |
| Methylphenidate (MPD) Test                   | Device | 50 ng/mL                         | /             |
| 6-Monoacetylmorphine (6-MAM) Test            | Device | 25/15/10/5/4 ng/mL               | CE            |
| Morphine (MOP) Test                          | Device | 15 ng/mL                         | CE            |
| Opiate (OPI) Test                            | Device | 50/40 ng/mL                      | CE            |
| Oxycodone (OXY) Test                         | Device | 50/40/20 ng/mL                   | CE            |
| Phencyclidine (PCP) Test                     | Device | 10 ng/mL                         | CE            |
| Phenytoin (PHEN) Test                        | Device | 150/100 ng/mL                    | /             |
| Pinaca Ab (K3) Test                          | Device | 10 ng/mL                         | /             |
| Pregabalin (PGB) Test                        | Device | 100 ng/mL                        | /             |
| Propoxyphene (PPX) Test                      | Device | 50/20 ng/mL                      | CE            |
| Synthetic Marijuana (K2) Test                | Device | 25/10/5 ng/mL                    | CE            |
| Tramadol (TRA) Test                          | Device | 100/50 ng/mL                     | CE            |
| Tricyclic Antidepressants (TCA) Test         | Device | 100 ng/mL                        | CE            |
| XLR-11 Test                                  | Device | 100 ng/mL                        | /             |
| Zolpidem (ZOL) Test                          | Device | 25 ng/mL                         | /             |
| Zopiclone (ZOP) Test                         | Device | 25 ng/mL                         | /             |
| Alcohol (ALC) Test                           | Device | 0.05/0.02%                       | CE            |

# NEOTEST™

## DOA SPLIT CUP URINE



### INSTRUCTIONS OF USE




# NEOTEST™

DOA URINE CUP



## INSTRUCTIONS OF USE




# NEOTEST™

DOA DIP CARD



## INSTRUCTIONS OF USE




# NEOTEST™

## DOA MULTI CASSETTE



## INSTRUCTIONS OF USE



## URINE DRUGS OF ABUSE (DOA) COMBINATIONS AVAILABLE

| Product Description                        | Format                      | Cut-off Value                           | Qualification |
|--------------------------------------------|-----------------------------|-----------------------------------------|---------------|
| Acetaminophen (ACE) Test                   | Strip/Cassette/Dip Card/Cup | 5000 ng/mL                              | CE            |
| Amphetamine (AMP) Test                     | Strip/Cassette/Dip Card/Cup | 2000/1000/500/300/250 ng/mL             | CE 510(k)     |
| Barbiturates (BAR) Test                    | Strip/Cassette/Dip Card/Cup | 2000/600/300/200 ng/mL                  | CE 510(k)     |
| Benzodiazepines (BZO) Test                 | Strip/Cassette/Dip Card/Cup | 600/400/300/200/100 ng/mL               | CE 510(k)     |
| Buprenorphine (BUP) Test                   | Strip/Cassette/Dip Card/Cup | 10/5 ng/mL                              | CE 510(k)     |
| Carisoprodol (SOMA) Test                   | Strip/Cassette/Dip Card/Cup | 1000 ng/mL                              | CE            |
| Clonazepam (CLO) Test                      | Strip/Cassette/Dip Card/Cup | 500/100 ng/mL                           | CE            |
| Cocaine (COC) Test                         | Strip/Cassette/Dip Card/Cup | 600/300/150/100 ng/mL                   | CE 510(k)     |
| Codeine (COD) Test                         | Strip/Cassette/Dip Card/Cup | 2000 ng/mL                              | CE            |
| Cotinine (COT) Test                        | Strip/Cassette/Dip Card/Cup | 400/300/200/100/50 ng/mL                | CE            |
| Ecstasy (MDMA) Test                        | Strip/Cassette/Dip Card/Cup | 2000/1000/500/300/250/150 ng/mL         | CE 510(k)     |
| Ethyl Glucuronide (EtG) Test               | Strip/Cassette/Dip Card/Cup | 500/300ng/mL                            | CE            |
| Fentanyl (FEN) Test                        | Strip/Cassette/Dip Card/Cup | 300/200/100/50 ng/mL                    | CE            |
| Norfentanyl (FEN) Test                     | Strip/Cassette/Dip Card/Cup | 200/50/20/10/5 ng/mL                    | CE            |
| Gabapentin (GAB) Test                      | Strip/Cassette/Dip Card/Cup | 3750/2000/1000 ng/mL                    | CE            |
| Hydrocodone (HCD) Test                     | Strip/Cassette/Dip Card/Cup | 300/10 ng/mL                            | CE            |
| Hydromorphone (HMO) Test                   | Strip/Cassette/Dip Card/Cup | 300 ng/mL                               | CE            |
| Ketamine (KET) Test                        | Strip/Cassette/Dip Card/Cup | 3000/2000/1000/500/100 ng/mL            | CE            |
| Kratom (KRA) Test                          | Strip/Cassette/Dip Card/Cup | 250/150/100 ng/mL                       | CE            |
| Lysergic acid diethylamide (LSD) Test      | Strip/Cassette/Dip Card/Cup | 20 ng/mL                                | CE            |
| Marijuana (THC) Test                       | Strip/Cassette/Dip Card/Cup | 600/300/200/150/50/40/25/20/18/15 ng/mL | CE 510(k)     |
| Methadone Metabolite (EDDP) Test           | Strip/Cassette/Dip Card/Cup | 300/100 ng/mL                           | CE 510(k)     |
| Methadone (MTD) Test                       | Strip/Cassette/Dip Card/Cup | 1000/600/300/200/50 ng/mL               | CE 510(k)     |
| Methamphetamine (MET) Test                 | Strip/Cassette/Dip Card/Cup | 2000/1000/500/300/250 ng/mL             | CE 510(k)     |
| Methaqualone (MQL) Test                    | Strip/Cassette/Dip Card/Cup | 300/1000 ng/mL                          | CE            |
| Methcathinone (MTC) Test                   | Strip/Cassette/Dip Card/Cup | 500/300 ng/mL                           | CE            |
| 3,4-Methylenedioxypyrovalerone (MDPV) Test | Strip/Cassette/Dip Card/Cup | 1000/500/300 ng/mL                      | CE            |
| Methylphenidate (MPD) Test                 | Strip/Cassette/Dip Card/Cup | 300 ng/mL                               | CE            |
| 6-Monoacetylmorphine (6-MAM) Test          | Strip/Cassette/Dip Card/Cup | 20/10 ng/mL                             | CE            |
| Morphine (MOP) Test                        | Strip/Cassette/Dip Card/Cup | 2000/600/300/150/100 ng/mL              | CE 510(k)     |
| Opiate (OPI) Test                          | Strip/Cassette/Dip Card/Cup | 2000/300/100 ng/mL                      | CE 510(k)     |
| Oxycodone (OXY) Test                       | Strip/Cassette/Dip Card/Cup | 300/100 ng/mL                           | CE 510(k)     |
| Phencyclidine (PCP) Test                   | Strip/Cassette/Dip Card/Cup | 50/25 ng/mL                             | CE 510(k)     |
| Pinaca Ab (K3) Test                        | Strip/Cassette/Dip Card/Cup | 10 ng/mL                                | CE            |
| Pregabalin (PGB) Test                      | Strip/Cassette/Dip Card/Cup | 2000/1000/500 ng/mL                     | CE            |
| Propoxyphene (PPX) Test                    | Strip/Cassette/Dip Card/Cup | 600/300 ng/mL                           | CE 510(k)     |
| Synthetic Marijuana (K2) Test              | Strip/Cassette/Dip Card/Cup | 75/50/25/20/10 ng/mL                    | CE            |
| Tramadol (TRA) Test                        | Strip/Cassette/Dip Card/Cup | 200/100 ng/mL                           | CE            |
| Tricyclic Antidepressants (TCA) Test       | Strip/Cassette/Dip Card/Cup | 1000/300 ng/mL                          | CE 510(k)     |
| UR-144 Test                                | Strip/Cassette/Dip Card/Cup | 50 ng/mL                                | CE            |
| Zolpidem (ZOL) Test                        | Strip/Cassette/Dip Card/Cup | 50 ng/mL                                | /             |
| Zopiclone (ZOP) Test                       | Strip/Cassette/Dip Card/Cup | 50 ng/mL                                | /             |
| Alcohol (ALC) Test                         | Strip/Cassette/Dip Card/Cup | 0.04%                                   | CE            |

| No | Device Name              | Model   | Product Identifier | Description                         |
|----|--------------------------|---------|--------------------|-------------------------------------|
| 1  | Drug of Abuse Rapid Test | NEOTEST | NT03020508         | Strips; Cocaine                     |
| 2  | Drug of Abuse Rapid Test | NEOTEST | NT03020608         | Strips; Amphetamine                 |
| 3  | Drug of Abuse Rapid Test | NEOTEST | NT03020808         | Strips; Methamphetamine             |
| 4  | Drug of Abuse Rapid Test | NEOTEST | NT03020708         | Strips; Barbiturates                |
| 5  | Drug of Abuse Rapid Test | NEOTEST | NT03021408         | Strips; Marijuana (K2)              |
| 6  | Drug of Abuse Rapid Test | NEOTEST | NT03021308         | Strips; Ketamine                    |
| 7  | Drug of Abuse Rapid Test | NEOTEST | NT03020908         | Strips; MDMA(Ecstasy)               |
| 8  | Drug of Abuse Rapid Test | NEOTEST | NT03021008         | Strips; Methadone                   |
| 9  | Drug of Abuse Rapid Test | NEOTEST | NT03021608         | Strips; Buprenorphine               |
| 10 | Drug of Abuse Rapid Test | NEOTEST | NT03020208         | Strips; Opiates                     |
| 11 | Drug of Abuse Rapid Test | NEOTEST | NT03021208         | Strips; PCP(Phencyclidine)          |
| 12 | Drug of Abuse Rapid Test | NEOTEST | NT03020208         | Strips; Tricyclic Antidepressants   |
| 13 | Drug of Abuse Rapid Test | NEOTEST | NT03020108         | Strips; Morphine                    |
| 14 | Drug of Abuse Rapid Test | NEOTEST | NT03021108         | Strips; Benzodiazaphine             |
| 15 | Drug of Abuse Rapid Test | NEOTEST | NT03020308         | Strips; Marijuana                   |
| 16 | Drug of Abuse Rapid Test | NEOTEST | NT03022308         | Strips; Alcohol (Urine)             |
| 17 | Drug of Abuse Rapid Test | NEOTEST | NT03022408         | Strips; KRATOM                      |
| 18 | Drug of Abuse Rapid Test | NEOTEST | NT03040508         | Cassette; Cocaine                   |
| 19 | Drug of Abuse Rapid Test | NEOTEST | NT03040608         | Cassette; Amphetamine               |
| 20 | Drug of Abuse Rapid Test | NEOTEST | NT03040808         | Cassette; Methamphetamine           |
| 21 | Drug of Abuse Rapid Test | NEOTEST | NT03040708         | Cassette; Barbiturates              |
| 22 | Drug of Abuse Rapid Test | NEOTEST | NT03041408         | Cassette; Marijuana (K2)            |
| 23 | Drug of Abuse Rapid Test | NEOTEST | NT03041308         | Cassette; Ketamine                  |
| 24 | Drug of Abuse Rapid Test | NEOTEST | NT03040908         | Cassette; MDMA(Ecstasy)             |
| 25 | Drug of Abuse Rapid Test | NEOTEST | NT03041008         | Cassette; Methadone                 |
| 26 | Drug of Abuse Rapid Test | NEOTEST | NT03041608         | Cassette; Buprenorphine             |
| 27 | Drug of Abuse Rapid Test | NEOTEST | NT03040208         | Cassette; Opiates                   |
| 28 | Drug of Abuse Rapid Test | NEOTEST | NT03041208         | Cassette; PCP(Phencyclidine)        |
| 29 | Drug of Abuse Rapid Test | NEOTEST | NT03040408         | Cassette; Tricyclic Antidepressants |
| 30 | Drug of Abuse Rapid Test | NEOTEST | NT03040108         | Cassette; Morphine                  |
| 31 | Drug of Abuse Rapid Test | NEOTEST | NT03041108         | Cassette; Benzodiazaphine           |
| 32 | Drug of Abuse Rapid Test | NEOTEST | NT03040308         | Cassette; Marijuana                 |
| 33 | Drug of Abuse Rapid Test | NEOTEST | NT03042308         | Cassette; Alcohol (Urine)           |
| 34 | Drug of Abuse Rapid Test | NEOTEST | NT03042408         | Cassette; KRATOM                    |
| 35 | Drug of Abuse Rapid Test | NEOTEST | NT0303040108       | Dip Card; COC MET OPI THC           |
| 36 | Drug of Abuse Rapid Test | NEOTEST | NT0303040208       | Dip Card; MOPTHCMET AMP             |
| 37 | Drug of Abuse Rapid Test | NEOTEST | NT0303050108       | Dip Card; AMP COC OPI PCP THC       |
| 38 | Drug of Abuse Rapid Test | NEOTEST | NT0303050208       | Dip Card; AMP COC MET OPI THC       |
| 39 | Drug of Abuse Rapid Test | NEOTEST | NT0303050308       | Dip Card; AMP COC MET MOP THC       |
| 40 | Drug of Abuse Rapid Test | NEOTEST | NT0303050408       | Dip Card; AMP COC MET MOP THC       |
| 41 | Drug of Abuse Rapid Test | NEOTEST | NT0303050508       | Dip Card; AMP MET MDMA OPI THC      |
| 42 | Drug of Abuse Rapid Test | NEOTEST | NT0303050608       | Dip Card; BZO COC MET OPI THC       |
| 43 | Drug of Abuse Rapid Test | NEOTEST | NT0303050708       | Dip Card; AMP BZO COC OPI THC       |
| 44 | Drug of Abuse Rapid Test | NEOTEST | NT0303050808       | Dip Card; MOP THC MET AMP KET       |
| 45 | Drug of Abuse Rapid Test | NEOTEST | NT0303050908       | Dip Card; AMP COC OPI KRA THC       |
| 46 | Drug of Abuse Rapid Test | NEOTEST | NT0303050908       | Dip Card; OPITHCCOCAMPK2            |
| 47 | Drug of Abuse Rapid Test | NEOTEST | NT0303051008       | Dip Card; MOPTHCCOCAMPK2            |
| 48 | Drug of Abuse Rapid Test | NEOTEST | NT0303060108       | Dip Card; AMP COC MET OPI PCP THC   |
| 49 | Drug of Abuse Rapid Test | NEOTEST | NT0303060208       | Dip Card; AMP BZO COC MET OPI THC   |
| 50 | Drug of Abuse Rapid Test | NEOTEST | NT0303060308       | Dip Card; AMP BZO COC MET MOP THC   |

| No | Device Name              | Model   | Product Identifier | Description                                                |
|----|--------------------------|---------|--------------------|------------------------------------------------------------|
| 51 | Drug of Abuse Rapid Test | NEOTEST | NT0303060408       | Dip Card; BZO COC MET OPI OXY THC                          |
| 52 | Drug of Abuse Rapid Test | NEOTEST | NT0303060508       | Dip Card; BZO COC MET MOP OXY THC                          |
| 53 | Drug of Abuse Rapid Test | NEOTEST | NT0303060608       | Dip Card; AMP COC MET OPI KRA THC                          |
| 54 | Drug of Abuse Rapid Test | NEOTEST | NT0303060708       | Dip Card; MOPTHCCOC mETBZOKET                              |
| 55 | Drug of Abuse Rapid Test | NEOTEST | NT0303060808       | Dip Card; MOPTHcCOCMETBZO K2                               |
| 56 | Drug of Abuse Rapid Test | NEOTEST | NT0303080108       | Dip Card; AMP BUP BZO COC MET OPI OXY THC                  |
| 57 | Drug of Abuse Rapid Test | NEOTEST | NT0303080208       | Dip Card; AMP BUP BZO COC MET OPI PCP THC                  |
| 58 | Drug of Abuse Rapid Test | NEOTEST | NT0303080308       | Dip Card; AMP BUP BZO COC MOP MET OXY THC                  |
| 59 | Drug of Abuse Rapid Test | NEOTEST | NT0303080408       | Dip Card; AMP BUP BZO COC MET OPI KRA THC                  |
| 60 | Drug of Abuse Rapid Test | NEOTEST | NT0303080508       | Dip Card; MOPTHCCOCAMPMETBZOKETOXY                         |
| 61 | Drug of Abuse Rapid Test | NEOTEST | NT0303080608       | Dip Card; MOPTHCCOCAMPMETBZOK2OXY                          |
| 62 | Drug of Abuse Rapid Test | NEOTEST | NT0303090108       | Dip Card; AMP BAR BZO COC MET MTD OPI PCP THC              |
| 63 | Drug of Abuse Rapid Test | NEOTEST | NT0303090208       | Dip Card; AMP BUP BZO COC MET MOP MDMA MTD THC             |
| 64 | Drug of Abuse Rapid Test | NEOTEST | NT0303100108       | Dip Card; AMP BAR BZO COC MET MTD OPI PCP TCA THC          |
| 65 | Drug of Abuse Rapid Test | NEOTEST | NT0303100208       | Dip Card; AMP BAR BUP BZO COC MET MTD OPI OXY THC          |
| 66 | Drug of Abuse Rapid Test | NEOTEST | NT0303100308       | Dip Card; AMP BAR BZO COC MET MOP MTD PCP TCA THC          |
| 67 | Drug of Abuse Rapid Test | NEOTEST | NT0303100408       | Dip Card; AMP BAR BUP BZO COC MET MOP MTD OXY THC          |
| 68 | Drug of Abuse Rapid Test | NEOTEST | NT0303100508       | Dip Card; AMP BAR BZO COC MET OPI MTD OXY PCP THC          |
| 69 | Drug of Abuse Rapid Test | NEOTEST | NT0303100608       | Dip Card; AMP BAR BZO COC MET MOP MTD OXY PCP THC          |
| 70 | Drug of Abuse Rapid Test | NEOTEST | NT0303100708       | Dip Card; BAR BUP BZO COC MET MTD OPI OXY TCA THC          |
| 71 | Drug of Abuse Rapid Test | NEOTEST | NT0303100808       | Dip Card; AMP BUP BZO COC MET MDMA MTD OPI OXY THC         |
| 72 | Drug of Abuse Rapid Test | NEOTEST | NT0303100908       | Dip Card; AMP BAR BZO COC MET OPI MTD OXY KRA THC          |
| 73 | Drug of Abuse Rapid Test | NEOTEST | NT0303101008       | Dip Card; OPITHCCOCAMPBUPMETMDMAMTDBZOKET                  |
| 74 | Drug of Abuse Rapid Test | NEOTEST | NT0303101108       | Dip Card; OPITHCCOCAMPBUPMETMDMAMTDBZK2                    |
| 75 | Drug of Abuse Rapid Test | NEOTEST | NT0303110108       | Dip Card; AMP BAR BZO COC MET MDMA MTD OPI OXY PCP THC     |
| 76 | Drug of Abuse Rapid Test | NEOTEST | NT0303110208       | Dip Card; AMP BAR BZO COC MET MDMA MOP MTD OXY PCPTHC      |
| 77 | Drug of Abuse Rapid Test | NEOTEST | NT0303120108       | Dip Card; AMP BAR BUP BZO COC MET MDMA MOP MTD OXY PCP THC |

| No  | Device Name              | Model   | Product Identifier | Description                                                |
|-----|--------------------------|---------|--------------------|------------------------------------------------------------|
| 78  | Drug of Abuse Rapid Test | NEOTEST | NT0303120208       | Dip Card; AMP BAR BUP BZO COC MDMA MTD TCA OPI OXY PCP THC |
| 79  | Drug of Abuse Rapid Test | NEOTEST | NT0303120308       | Dip Card; AMP BAR BZO COC MET MDMA MOP MTD OXY PPX TCA THC |
| 80  | Drug of Abuse Rapid Test | NEOTEST | NT0303120408       | Dip Card; AMP BAR BUP BZO COC MET MDMA MTD OPI OXY PCP THC |
| 81  | Drug of Abuse Rapid Test | NEOTEST | NT0303120508       | Dip Card; OPITHCCOCAMPBARMETMDMAMTDBZOK ETK2OXY            |
| 82  | Drug of Abuse Rapid Test | NEOTEST | NT0301050108       | Cup; AMP COC OPI PCP THC                                   |
| 83  | Drug of Abuse Rapid Test | NEOTEST | NT0301050208       | Cup; AMP COC MET OPI THC                                   |
| 84  | Drug of Abuse Rapid Test | NEOTEST | NT0301050308       | Cup; AMP COC MET MOP THC                                   |
| 85  | Drug of Abuse Rapid Test | NEOTEST | NT0301050408       | Cup; AMP COC MET MOP THC                                   |
| 86  | Drug of Abuse Rapid Test | NEOTEST | NT0301050508       | Cup; AMP MET MDMA OPI THC                                  |
| 87  | Drug of Abuse Rapid Test | NEOTEST | NT0301050608       | Cup; BZO COC MET OPI THC                                   |
| 88  | Drug of Abuse Rapid Test | NEOTEST | NT0301050708       | Cup; AMP BZO COC OPI THC                                   |
| 89  | Drug of Abuse Rapid Test | NEOTEST | NT0301050808       | Cup; MOP THC MET AMP KET                                   |
| 90  | Drug of Abuse Rapid Test | NEOTEST | NT0301050908       | Cup; AMP COC OPI KRA THC                                   |
| 91  | Drug of Abuse Rapid Test | NEOTEST | NT0301051008       | Cup; OPITHCCOCAMPK2                                        |
| 92  | Drug of Abuse Rapid Test | NEOTEST | NT0301051108       | Cup; MOPTHCCOCAMPK2                                        |
| 93  | Drug of Abuse Rapid Test | NEOTEST | NT0301060108       | Cup; AMP COC MET OPI PCP THC                               |
| 94  | Drug of Abuse Rapid Test | NEOTEST | NT0301060208       | Cup;AMP BZO COC MET OPI THC                                |
| 95  | Drug of Abuse Rapid Test | NEOTEST | NT0301060308       | Cup;AMP BZO COCMET MOP THC                                 |
| 96  | Drug of Abuse Rapid Test | NEOTEST | NT0301060408       | Cup; BZO COC MET OPI KET THC                               |
| 97  | Drug of Abuse Rapid Test | NEOTEST | NT0301060508       | Cup;BZO COC MET MOP OXY THC                                |
| 98  | Drug of Abuse Rapid Test | NEOTEST | NT0301060608       | Cup; AMP COC MET OPI KRA THC                               |
| 99  | Drug of Abuse Rapid Test | NEOTEST | NT0301060708       | Cup;OPITHCCOCAMP MDMA BZO                                  |
| 100 | Drug of Abuse Rapid Test | NEOTEST | NT0301060808       | Cup;AMP COC OPI PCP THCALC                                 |
| 101 | Drug of Abuse Rapid Test | NEOTEST | NT0301060908       | Cup;AMP COC MET OPI THC ALC                                |
| 102 | Drug of Abuse Rapid Test | NEOTEST | NT0301061008       | Cup;AMP COC MET MOP THC ALC                                |
| 103 | Drug of Abuse Rapid Test | NEOTEST | NT0301061108       | Cup;AMP COC MET MOP THC ALC                                |
| 104 | Drug of Abuse Rapid Test | NEOTEST | NT0301061208       | Cup;AMP MET MDMA OPI THC ALC                               |
| 105 | Drug of Abuse Rapid Test | NEOTEST | NT0301061308       | Cup;BZO COC MET OPI THC ALC                                |
| 106 | Drug of Abuse Rapid Test | NEOTEST | NT0301061408       | Cup;AMP BZO COC OPI THCALC                                 |
| 107 | Drug of Abuse Rapid Test | NEOTEST | NT0301061508       | Cup;MOP THC MET AMP KETALC                                 |
| 108 | Drug of Abuse Rapid Test | NEOTEST | NT0301061608       | Cup;AMP COC OPI KRA THC ALC                                |
| 109 | Drug of Abuse Rapid Test | NEOTEST | NT0301061708       | Cup;AMP COC OPI PCP THCADL                                 |
| 110 | Drug of Abuse Rapid Test | NEOTEST | NT0301061808       | Cup;AMP COC MET OPI THC ADL                                |
| 111 | Drug of Abuse Rapid Test | NEOTEST | NT0301061908       | Cup;AMP COC MET MOP THC ADL                                |
| 112 | Drug of Abuse Rapid Test | NEOTEST | NT0301062008       | Cup;AMP COC MET MOP THC ADL                                |
| 113 | Drug of Abuse Rapid Test | NEOTEST | NT0301062108       | Cup;AMP MET MDMA OPI THC ADL                               |
| 114 | Drug of Abuse Rapid Test | NEOTEST | NT0301062208       | Cup;BZO COC MET OPI THC ADL                                |
| 115 | Drug of Abuse Rapid Test | NEOTEST | NT0301062308       | Cup;AMP BZO COC OPI THCADL                                 |
| 116 | Drug of Abuse Rapid Test | NEOTEST | NT0301062408       | Cup;MOP THC MET AMP KETADL                                 |
| 117 | Drug of Abuse Rapid Test | NEOTEST | NT0301062508       | Cup;AMP COC OPI KRA THC ADL                                |
| 118 | Drug of Abuse Rapid Test | NEOTEST | NT0301062608       | Cup;OPITHCCOCAMPK2ADL                                      |
| 119 | Drug of Abuse Rapid Test | NEOTEST | NT0301062708       | Cup;MOPTHCCOCAMPK2ADL                                      |
| 120 | Drug of Abuse Rapid Test | NEOTEST | NT0301062808       | Cup;OPITHCCOCAMPK2ACL                                      |
| 121 | Drug of Abuse Rapid Test | NEOTEST | NT0301062908       | Cup;MOPTHCCOCAMPK2ACL                                      |
| 122 | Drug of Abuse Rapid Test | NEOTEST | NT0301070108       | Cup;OPI THC COCAMPMET MDMA BZO                             |
| 123 | Drug of Abuse Rapid Test | NEOTEST | NT0301070208       | Cup;AMP COC OPI PCP THCALCADL                              |
| 124 | Drug of Abuse Rapid Test | NEOTEST | NT0301070308       | Cup;AMP COC MET OPI THC ALCADL                             |

| No  | Device Name              | Model   | Product Identifier | Description                              |
|-----|--------------------------|---------|--------------------|------------------------------------------|
| 125 | Drug of Abuse Rapid Test | NEOTEST | NT0301070408       | Cup;AMP COC MET MOP THC ALCADL           |
| 126 | Drug of Abuse Rapid Test | NEOTEST | NT0301070508       | Cup;AMP COC MET MOP THC ALCADL           |
| 127 | Drug of Abuse Rapid Test | NEOTEST | NT0301070608       | Cup;AMP MET MDMA OPI THC ALCADL          |
| 128 | Drug of Abuse Rapid Test | NEOTEST | NT0301070708       | Cup;BZO COC MET OPI THC ALCADL           |
| 129 | Drug of Abuse Rapid Test | NEOTEST | NT0301070808       | Cup;AMP BZO COC OPI THCALCADL            |
| 130 | Drug of Abuse Rapid Test | NEOTEST | NT0301070908       | Cup;MOP THC MET AMP KETALCADL            |
| 131 | Drug of Abuse Rapid Test | NEOTEST | NT0301071008       | Cup;AMP COC OPI KRA THC ALCADL           |
| 132 | Drug of Abuse Rapid Test | NEOTEST | NT0301071108       | Cup;OPITHCCOCAMPK2ADLACL                 |
| 133 | Drug of Abuse Rapid Test | NEOTEST | NT0301071208       | Cup;MOPTHCCOCAMPK2ADLACL                 |
| 134 | Drug of Abuse Rapid Test | NEOTEST | NT0301080108       | Cup;AMP BUP BZO COC MET OPI OXY THC      |
| 135 | Drug of Abuse Rapid Test | NEOTEST | NT0301080208       | Cup;AMP BUP BZO COC MET OPI PCP THC      |
| 136 | Drug of Abuse Rapid Test | NEOTEST | NT0301080308       | Cup;AMP BUP BZO COC MOP MET OXY THC      |
| 137 | Drug of Abuse Rapid Test | NEOTEST | NT0301080408       | Cup;AMP BUP BZO COC MET OPI KRA THC      |
| 138 | Drug of Abuse Rapid Test | NEOTEST | NT0301080508       | Cup;OPI THC COCAMPMET MDMA BZOALC        |
| 139 | Drug of Abuse Rapid Test | NEOTEST | NT0301080608       | Cup;MOPTHCCOCAMPMETMDMABZOKET            |
| 140 | Drug of Abuse Rapid Test | NEOTEST | NT0301080708       | Cup;MOPTHCCOCAMPMETMDMABZOK2             |
| 141 | Drug of Abuse Rapid Test | NEOTEST | NT0301090108       | Cup;AMP BAR BZO COC MET MTD OPI PCP THC  |
| 142 | Drug of Abuse Rapid Test | NEOTEST | NT0301090208       | Cup;AMP BUP BZO COC MET MOP MDMA MTD THC |
| 143 | Drug of Abuse Rapid Test | NEOTEST | NT0301090308       | Cup;OPI THC COCAMPMET MDMA BZOALCADL     |
| 144 | Drug of Abuse Rapid Test | NEOTEST | NT0301090408       | Cup;AMP BUP BZO COC MET OPI OXY THC ALC  |
| 145 | Drug of Abuse Rapid Test | NEOTEST | NT0301090508       | Cup;AMP BUP BZO COC MET OPI PCP THC ALC  |
| 146 | Drug of Abuse Rapid Test | NEOTEST | NT0301090608       | Cup;AMP BUP BZO COC MOP MET OXY THCALC   |
| 147 | Drug of Abuse Rapid Test | NEOTEST | NT0301090708       | Cup;AMP BUP BZO COC MET OPI KRA THC ALC  |
| 148 | Drug of Abuse Rapid Test | NEOTEST | NT0301090808       | Cup;AMP BUP BZO COC MET OPI OXY THC ADL  |
| 149 | Drug of Abuse Rapid Test | NEOTEST | NT0301090908       | Cup;AMP BUP BZO COC MET OPI PCP THC ADL  |
| 150 | Drug of Abuse Rapid Test | NEOTEST | NT0301091008       | Cup;AMP BUP BZO COC MOP MET OXY THCADL   |
| 151 | Drug of Abuse Rapid Test | NEOTEST | NT0301091108       | Cup;AMP BUP BZO COC MET OPI KRA THC ADL  |
| 152 | Drug of Abuse Rapid Test | NEOTEST | NT0301091208       | Cup;MOPTHCCOCAMPMETMDMABZOKETA DL        |
| 153 | Drug of Abuse Rapid Test | NEOTEST | NT0301091308       | Cup;MOPTHCCOCAMPMETMDMABZOK2A DL         |
| 154 | Drug of Abuse Rapid Test | NEOTEST | NT0301091408       | Cup;MOPTHCCOCAMPMETMDMABZOKETA LC        |

| No  | Device Name              | Model   | Product Identifier | Description                                       |
|-----|--------------------------|---------|--------------------|---------------------------------------------------|
| 155 | Drug of Abuse Rapid Test | NEOTEST | NT0301091508       | Cup;MOPTHCCOCAMPMETMDMABZOK2ALC                   |
| 156 | Drug of Abuse Rapid Test | NEOTEST | NT0301100108       | Cup;AMP BAR BZO COC MET MTD OPI PCP TCA THC       |
| 157 | Drug of Abuse Rapid Test | NEOTEST | NT0301100208       | Cup; AMP BAR BUP BZO COC MET MTD OPI OXY THC      |
| 158 | Drug of Abuse Rapid Test | NEOTEST | NT0301100308       | Cup;AMP BAR BZO COC MET MOP MTD PCP TCA THC       |
| 159 | Drug of Abuse Rapid Test | NEOTEST | NT0301100408       | Cup;AMP BAR BUP BZO COC MET MOP MTD OXY THC       |
| 160 | Drug of Abuse Rapid Test | NEOTEST | NT0301100508       | Cup;AMP BAR BZO COC MET OPI MTD OXY PCP THC       |
| 161 | Drug of Abuse Rapid Test | NEOTEST | NT0301100608       | Cup;AMP BAR BZO COC MET MOP MTD OXY PCP THC       |
| 162 | Drug of Abuse Rapid Test | NEOTEST | NT0301100708       | Cup;BAR BUP BZO COC MET MTD OPI OXY TCA THC       |
| 163 | Drug of Abuse Rapid Test | NEOTEST | NT0301100808       | Cup;AMP BUP BZO COC MET MDMA MTD OPI OXY THC      |
| 164 | Drug of Abuse Rapid Test | NEOTEST | NT0301100908       | Cup;AMP BAR BZO COC MET OPI MTD OXY KRA THC       |
| 165 | Drug of Abuse Rapid Test | NEOTEST | NT0301101008       | Cup;AMP BAR BZO COC MET MTD OPI PCP THCALC        |
| 166 | Drug of Abuse Rapid Test | NEOTEST | NT0301101108       | Cup;AMP BUP BZO COC MET MOP MDMA MTD THCADL       |
| 167 | Drug of Abuse Rapid Test | NEOTEST | NT0301101108       | Cup;MOPTHCCOCAMPMETMDMAMTDBZOK2BUP                |
| 168 | Drug of Abuse Rapid Test | NEOTEST | NT0301101208       | Cup;OPITHCCOCAMPMETMDMAMTDBZOKETBUP               |
| 169 | Drug of Abuse Rapid Test | NEOTEST | NT0301110108       | Cup; AMP BAR BZO COC MET MDMA MTD OPI OXY PCP THC |
| 170 | Drug of Abuse Rapid Test | NEOTEST | NT0301110208       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PCPTHC   |
| 171 | Drug of Abuse Rapid Test | NEOTEST | NT0301110308       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY TCA THC  |
| 172 | Drug of Abuse Rapid Test | NEOTEST | NT0301110408       | Cup;AMP BAR BZO COC MET MTD OPI PCP THCALCADL     |
| 173 | Drug of Abuse Rapid Test | NEOTEST | NT0301110508       | Cup;AMP BAR BZO COC MET MTD OPI PCP TCA THC ALC   |
| 174 | Drug of Abuse Rapid Test | NEOTEST | NT0301110608       | Cup; AMP BAR BUP BZO COC MET MTD OPI OXY THC ALC  |
| 175 | Drug of Abuse Rapid Test | NEOTEST | NT0301110708       | Cup;AMP BAR BZO COC MET MOP MTD PCP TCA THC ALC   |
| 176 | Drug of Abuse Rapid Test | NEOTEST | NT0301110808       | Cup;AMP BAR BUP BZO COC MET MOP MTD OXY THC ALC   |
| 177 | Drug of Abuse Rapid Test | NEOTEST | NT0301110908       | Cup;AMP BAR BZO COC MET OPI MTD OXY PCP THC ALC   |
| 178 | Drug of Abuse Rapid Test | NEOTEST | NT0301111008       | Cup;AMP BAR BZO COC MET MOP MTD OXY PCP THC ALC   |
| 179 | Drug of Abuse Rapid Test | NEOTEST | NT0301111108       | Cup;BAR BUP BZO COC MET MTD OPI OXY TCA THC ALC   |
| 180 | Drug of Abuse Rapid Test | NEOTEST | NT0301111208       | Cup;AMP BUP BZO COC MET MDMA MTD OPI OXY THCALC   |
| 181 | Drug of Abuse Rapid Test | NEOTEST | NT0301111308       | Cup;AMP BAR BZO COC MET OPI MTD OXY KRA THC ALC   |
| 182 | Drug of Abuse Rapid Test | NEOTEST | NT0301111408       | Cup;AMP BAR BZO COC MET MTD OPI PCP TCA THC ADL   |

| No  | Device Name              | Model   | Product Identifier | Description                                           |
|-----|--------------------------|---------|--------------------|-------------------------------------------------------|
| 183 | Drug of Abuse Rapid Test | NEOTEST | NT0301111508       | Cup; AMP BAR BUP BZO COC MET MTD OPI OXY THC ADL      |
| 184 | Drug of Abuse Rapid Test | NEOTEST | NT0301111608       | Cup;AMP BAR BZO COC MET MOP MTD PCP TCA THC ADL       |
| 185 | Drug of Abuse Rapid Test | NEOTEST | NT0301111708       | Cup;AMP BAR BUP BZO COC MET MOP MTD OXY THC ADL       |
| 186 | Drug of Abuse Rapid Test | NEOTEST | NT0301111808       | Cup;AMP BAR BZO COC MET OPI MTD OXY PCP THC ADL       |
| 187 | Drug of Abuse Rapid Test | NEOTEST | NT0301111908       | Cup;AMP BAR BZO COC MET MOP MTD OXY PCP THC ADL       |
| 188 | Drug of Abuse Rapid Test | NEOTEST | NT0301112008       | Cup;BAR BUP BZO COC MET MTD OPI OXY TCA THC ADL       |
| 189 | Drug of Abuse Rapid Test | NEOTEST | NT0301112108       | Cup;AMP BUP BZO COC MET MDMA MTD OPI OXY THCADL       |
| 190 | Drug of Abuse Rapid Test | NEOTEST | NT0301112208       | Cup;AMP BAR BZO COC MET OPI MTD OXY KRA THC ADL       |
| 191 | Drug of Abuse Rapid Test | NEOTEST | NT0301112308       | Cup;MOPTHCCOCAMPMETMDMAMTDBZO K2BUPADL                |
| 192 | Drug of Abuse Rapid Test | NEOTEST | NT0301112408       | Cup;OPITHCCOCAMPMETMDMAMTDBZOK ETBUPADL               |
| 193 | Drug of Abuse Rapid Test | NEOTEST | NT0301112508       | Cup;MOPTHCCOCAMPMETMDMAMTDBZO K2BUPACL                |
| 194 | Drug of Abuse Rapid Test | NEOTEST | NT0301112608       | Cup;OPITHCCOCAMPMETMDMAMTDBZOK ETBUPACL               |
| 195 | Drug of Abuse Rapid Test | NEOTEST | NT0301120108       | Cup;AMP BAR BUP BZO COC MET MDMA MOP MTD OXY PCP THC  |
| 196 | Drug of Abuse Rapid Test | NEOTEST | NT0301120208       | Cup;AMP BAR BUP BZO COC MDMA MTD TCA OPI OXY PCP THC  |
| 197 | Drug of Abuse Rapid Test | NEOTEST | NT0301120308       | Cup; AMP BAR BUP BZO COC MET MDMA MTD OPI OXY PCP THC |
| 198 | Drug of Abuse Rapid Test | NEOTEST | NT0301120408       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PPX TCA THC  |
| 199 | Drug of Abuse Rapid Test | NEOTEST | NT0301120508       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY KRA TCA THC  |
| 200 | Drug of Abuse Rapid Test | NEOTEST | NT0301120608       | Cup; AMP BAR BZO COC MET MDMA MTD OPI OXY PCP THC ALC |
| 201 | Drug of Abuse Rapid Test | NEOTEST | NT0301120708       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PCPTHC ALC   |
| 202 | Drug of Abuse Rapid Test | NEOTEST | NT0301120808       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY TCA THC ALC  |
| 203 | Drug of Abuse Rapid Test | NEOTEST | NT0301120908       | Cup; AMP BAR BZO COC MET MDMA MTD OPI OXY PCP THC ADL |
| 204 | Drug of Abuse Rapid Test | NEOTEST | NT0301121008       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PCPTHC ADL   |
| 205 | Drug of Abuse Rapid Test | NEOTEST | NT0301121108       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY TCA THC ADL  |
| 206 | Drug of Abuse Rapid Test | NEOTEST | NT0301121208       | Cup;AMP BAR BZO COC MET MTD OPI PCP TCA THC ALCADL    |
| 207 | Drug of Abuse Rapid Test | NEOTEST | NT0301121308       | Cup; AMP BAR BUP BZO COC MET MTD OPI OXY THC ALCADL   |
| 208 | Drug of Abuse Rapid Test | NEOTEST | NT0301121408       | Cup;AMP BAR BZO COC MET MOP MTD PCP TCA THC ALCADL    |
| 209 | Drug of Abuse Rapid Test | NEOTEST | NT0301121508       | Cup;AMP BAR BUP BZO COC MET MOP MTD OXY THC ALCADL    |
| 210 | Drug of Abuse Rapid Test | NEOTEST | NT0301121608       | Cup;AMP BAR BZO COC MET OPI MTD OXY PCP THC ALCADL    |

| No  | Device Name              | Model   | Product Identifier | Description                                                 |
|-----|--------------------------|---------|--------------------|-------------------------------------------------------------|
| 211 | Drug of Abuse Rapid Test | NEOTEST | NT0301121708       | Cup;AMP BAR BZO COC MET MOP MTD OXY PCP THC ALCADL          |
| 212 | Drug of Abuse Rapid Test | NEOTEST | NT0301121808       | Cup;BAR BUP BZO COC MET MTD OPI OXY TCA THC ALCADL          |
| 213 | Drug of Abuse Rapid Test | NEOTEST | NT0301121908       | Cup;AMP BUP BZO COC MET MDMA MTD OPI OXY THCALCADL          |
| 214 | Drug of Abuse Rapid Test | NEOTEST | NT0301122008       | Cup;AMP BAR BZO COC MET OPI MTD OXY KRA THC ALCADL          |
| 215 | Drug of Abuse Rapid Test | NEOTEST | NT0301122108       | Cup;MOPTHCCOCAMPMETMDMAMTDBZO K2BUPADLACL                   |
| 216 | Drug of Abuse Rapid Test | NEOTEST | NT0301122208       | Cup;OPITHCCOCAMPMETMDMAMTDBZOK ETBUPADLACL                  |
| 217 | Drug of Abuse Rapid Test | NEOTEST | NT0301122308       | Cup;OPITHCCOCAMPBARMETMDMAMTDB ZOKETK2OXY                   |
| 218 | Drug of Abuse Rapid Test | NEOTEST | NT0301130108       | Cup;AMP BAR BUP BZO COC MET MDMA MTD MOP OXY PCP TCA THC    |
| 219 | Drug of Abuse Rapid Test | NEOTEST | NT0301130208       | Cup;AMP BAR BUP BZO COC MET MDMA MOP MTD OXY PCP THC ALC    |
| 220 | Drug of Abuse Rapid Test | NEOTEST | NT0301130308       | Cup;AMP BAR BUP BZO COC MDMA MTD TCA OPI OXY PCP THC ALC    |
| 221 | Drug of Abuse Rapid Test | NEOTEST | NT0301130408       | Cup; AMP BAR BUP BZO COC MET MDMA MTD OPI OXY PCP THCALC    |
| 222 | Drug of Abuse Rapid Test | NEOTEST | NT0301130508       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PPX TCA THC ALC    |
| 223 | Drug of Abuse Rapid Test | NEOTEST | NT0301130608       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY KRA TCA THC ALC    |
| 224 | Drug of Abuse Rapid Test | NEOTEST | NT0301130708       | Cup;AMP BAR BUP BZO COC MET MDMA MOP MTD OXY PCP THC ADL    |
| 225 | Drug of Abuse Rapid Test | NEOTEST | NT0301130808       | Cup;AMP BAR BUP BZO COC MDMA MTD TCA OPI OXY PCP THC ADL    |
| 226 | Drug of Abuse Rapid Test | NEOTEST | NT0301130908       | Cup; AMP BAR BUP BZO COC MET MDMA MTD OPI OXY PCP THCADL    |
| 227 | Drug of Abuse Rapid Test | NEOTEST | NT0301131008       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY PPX TCA THC ADL    |
| 228 | Drug of Abuse Rapid Test | NEOTEST | NT0301131108       | Cup;AMP BAR BZO COC MET MDMA MOP MTD OXY KRA TCA THC ADL    |
| 229 | Drug of Abuse Rapid Test | NEOTEST | NT0301131208       | Cup;OPITHCCOCAMPBARMETMDMAMTDB ZOKETK2OXYADL                |
| 230 | Drug of Abuse Rapid Test | NEOTEST | NT0301131308       | Cup;OPITHCCOCAMPBARMETMDMAMTDB ZOKETK2OXYALC                |
| 231 | Drug of Abuse Rapid Test | NEOTEST | NT0301140108       | Cup;AMP BAR BUP BZO COC MET MDMA MTD MOP OXY PCP TCA THCALC |
| 232 | Drug of Abuse Rapid Test | NEOTEST | NT0301140208       | Cup;AMP BAR BUP BZO COC MET MDMA MTD MOP OXY PCP TCA THCADL |
| 233 | Drug of Abuse Rapid Test | NEOTEST | NT0301140308       | Cup;AMP BAR BUP BZO COC MET MDMA MOP MTD OXY PCP THC ALCADL |
| 234 | Drug of Abuse Rapid Test | NEOTEST | NT0301140408       | Cup;AMP BAR BUP BZO COC MDMA MTD TCA OPI OXY PCP THC ALCADL |
| 235 | Drug of Abuse Rapid Test | NEOTEST | NT0301140508       | Cup; AMP BAR BUP BZO COC MET MDMA MTD OPI OXY PCP THCALCADL |

| No  | Device Name              | Model   | Product Identifier | Description                                                    |
|-----|--------------------------|---------|--------------------|----------------------------------------------------------------|
| 236 | Drug of Abuse Rapid Test | NEOTEST | NT0301140608       | Cup;AMP BAR BZO COC MET MDMA MOP<br>MTD OXY PPX TCA THC ALCADL |
| 237 | Drug of Abuse Rapid Test | NEOTEST | NT0301140708       | Cup;AMP BAR BZO COC MET MDMA MOP<br>MTD OXY KRA TCA THC ALCADL |
| 238 | Drug of Abuse Rapid Test | NEOTEST | NT0301140808       | Cup;OPITHCCOCAMPBARMETMDMAMTDB<br>ZOKETK2OXYADLALC             |